

## TRIGEN<sup>◇</sup> INTERTAN<sup>◇</sup> Intertrochanteric Antegrade Nail: clinical and economic benefits in the treatment of intertrochanteric hip fractures

### Discussion points

- Treatment of intertrochanteric hip fractures represents a substantial proportion of the high overall clinical and economic burden of hip fracture<sup>1,2</sup>
- TRIGEN INTERTAN has shown improved biomechanical performance<sup>3-5</sup> and clinical outcomes, including revision/reoperation rate,<sup>6</sup> compared with other intramedullary (IM) nails in the treatment of intertrochanteric fractures
- In an economic analysis, TRIGEN INTERTAN has been associated with reduced per patient hospital costs, when compared with alternative IM nails<sup>7</sup>

### The clinical and economic burden of hip fractures

Though the incidence of hip fractures is relatively low, accounting for 15% of fractures treated, they are associated with high costs;<sup>8</sup> for example, in the Netherlands, hip fractures have been estimated to contribute 53% of osteoporosis-related fracture costs.<sup>8</sup> Hip fractures have a high mortality rate, even when surgically treated,<sup>9-11</sup> and are associated with high rates of morbidity, with considerable potential impact on patients' mobility and daily living.<sup>12,13</sup> The likelihood that a hip fracture patient will recover to their pre-fracture level of function is less than 50%.<sup>13</sup> Moreover, many patients require reoperation due to complications,<sup>10</sup> contributing to the already substantial economic burden associated with hip fractures (Figure 1).<sup>14,15</sup>




Figure 1. The clinical and economic burden of hip fracture

### Intertrochanteric hip fractures

Intertrochanteric fractures are hip fractures that occur between the greater and lesser trochanters of the proximal femur.<sup>14</sup> They comprise a large proportion of both the overall clinical burden and total direct medical costs relating to hip fractures (Figure 2).<sup>1</sup>




Figure 2. The contribution of intertrochanteric hip fracture to overall hip fracture burden

## + Evidence in focus

### Treatment of intertrochanteric hip fractures

The standard of care for intertrochanteric fractures is surgical fixation, unless contraindicated,<sup>16</sup> with the aim of achieving early mobilisation and facilitating earlier rehabilitation and functional recovery.<sup>17</sup> A number of different surgical treatment approaches can be used, with the procedure primarily influenced by fracture pattern.<sup>16,18</sup>


- Intramedullary nails are the most commonly used fixation method for intertrochanteric fractures in North America,<sup>19</sup> with a range of types available including TRIGEN<sup>®</sup> INTERTAN<sup>®</sup> Intertrochanteric Antegrade Nail, Gamma3<sup>™</sup> Nailing System (Stryker, Michigan, USA), TFN-ADVANCED<sup>™</sup> Proximal Femoral Nailing System (TFNA<sup>™</sup>; DePuy Synthes, Pennsylvania, USA) and Proximal Femoral Nail Antirotation (PFNA<sup>™</sup>; DePuy Synthes, Solothurn, Switzerland)
- Dynamic or sliding hip screws may be used for stable fractures but have a greater risk of failure in older patients with a lower bone density,<sup>18,20</sup> and have limitations in patients with unstable fractures<sup>21,22</sup>
- Hemiarthroplasty, an operation where the femoral head alone is replaced with a prosthesis, typically has higher complication rates but may be used for a small number of frail, elderly patients.<sup>18</sup> For example, it is recommended for patients with ipsilateral hip osteoarthritis, unstable fracture patterns where bone quality is poor, and failed internal fixations<sup>18</sup>

### TRIGEN INTERTAN: An IM nail offering improved performance over comparator devices

TRIGEN INTERTAN differs from traditional single-screw IM devices in its use of two Integrated Compression Screws and in other aspects of its design:

- The dual-screw design is intended to maintain the stability of the femoral head by pulling it tightly to the hip bone and preventing it from rotating out of position<sup>5,23</sup>
- The trapezoidal shape of the implant is designed to provide a tighter fit in the hip bone and therefore additional strength where forces tend to be highest<sup>24</sup>
- Finally, the lower tip of the nail is less rigid, and designed to reduce stress on the hip bone


TRIGEN INTERTAN has been shown to result in improved performance in biomechanical testing compared with other IM nails (Figure 3).<sup>3-5</sup>



### TRIGEN INTERTAN has the potential to improve clinical outcomes in unstable intertrochanteric fractures

A meta-analysis of studies reporting on adult patients with unstable intertrochanteric fractures has permitted the direct comparison of clinical and functional outcomes with TRIGEN INTERTAN and other IM nails including Gamma3<sup>™</sup> and PFNA<sup>™</sup> (Figure 4).<sup>6</sup>

TRIGEN INTERTAN demonstrated reductions in the risk of revision/reoperation, implant-related failure, and hip and thigh pain when compared with other IM nails (Figure 4). Together, these findings highlight the opportunity to improve outcomes for patients through the use of TRIGEN INTERTAN in unstable intertrochanteric fractures.



## + Evidence in focus

### TRIGEN<sup>◇</sup> INTERTAN<sup>◇</sup> reduces hospital costs compared with other IM nails for the treatment of intertrochanteric hip fractures



The US Premier database comprises data from over 700 hospitals with approximately one billion patient encounters,<sup>25</sup> and is broadly representative of the US hospital experience.<sup>26</sup> Using this database, all inpatient costs were compared for patients who:

- Were admitted with unstable intertrochanteric fractures – identified using ICD-10 diagnosis code S72.14 and a recorded Diagnosis-Related Group (DRG) of 480, 481 or 482 (Figure 5) – between January 2017 and September 2019
- Were treated using either TRIGEN INTERTAN (n=1,585), Gamma3™ (n=5,538) or TFNA™ nails (n=7,078)

Departmental costs were evaluated for the most common of the three DRGs (481), which includes patients with hip or femur fractures (except major joint) and complications or comorbidities. This group excludes patients who do not have complications or comorbidities, as well as those with major complications and comorbidities, and thus, overall, represents cases of moderate severity.<sup>27</sup>

In a matched comparison of Gamma3™ and TRIGEN INTERTAN nails, TRIGEN INTERTAN was associated with a mean \$1,431 reduction in hospital costs per patient (Figure 6a).<sup>7</sup> TRIGEN INTERTAN also resulted in mean hospital cost savings of \$968 per patient in a matched comparison with TFNA™ nail (Figure 6a).<sup>7</sup> Furthermore, this analysis highlighted that the cost of the IM nail implants themselves represent <15% of the total hospital costs (Figure 6b).<sup>7</sup>

TRIGEN INTERTAN was also cost saving compared with Gamma3™ and TFNA™ for DRGs 480 and 482. Mean hospital cost savings compared with Gamma3™ were \$1,637 and \$1,176 for DRGs 480 and 482, respectively. Similarly, mean hospital cost savings compared with TFNA™ were \$543 and \$693 for DRGs 480 and 482, respectively.<sup>7</sup>



### Conclusions

- Intertrochanteric fractures represent a **substantial proportion of the overall clinical and economic burden of hip fracture**<sup>1</sup>
- The design features of TRIGEN<sup>®</sup> INTERTAN<sup>®</sup> help to improve **biomechanical strength**<sup>3–5</sup> and **stability**<sup>3,5</sup> over other IM nails
- TRIGEN INTERTAN is associated with **significant improvements in clinical outcomes compared with other IM nails**, including revision/reoperation rate, implant failure and hip and thigh pain<sup>6</sup>
- TRIGEN INTERTAN has the potential to generate **substantial savings by reducing hospital costs**<sup>7</sup>

### References

1. Adeyemi A, Delhougne G. Incidence and economic burden of intertrochanteric fracture: a Medicare claims database analysis. *JB JS Open Access*. 2019;4:e0045.
2. Veronese N, Maggi S. Epidemiology and social costs of hip fracture. *Injury*. 2018;49:1458–1460.
3. Huang Y, Zhang C, Luo Y. A comparative biomechanical study of proximal femoral nail (InterTAN) and proximal femoral nail antirotation for intertrochanteric fractures. *Int Orthop*. 2013;37:2465–73.
4. Kubiak EN, Bong M, Park SS, et al. Intramedullary fixation of unstable intertrochanteric hip fractures: one or two lag screws. *J Orthop Trauma*. 2004;18:12–7.
5. Santoni BG, Nayak AN, Cooper SA, et al. Comparison of femoral head rotation and varus collapse between a single lag screw and integrated dual screw intertrochanteric hip fracture fixation device using a cadaveric hemi-pelvis biomechanical model. *J Orthop Trauma*. 2016;30:164–9.
6. Smith+Nephew 2021. EA\_TRAUMA\_INTERTAN\_003\_v1. 2021.
7. Smith+Nephew. Evidence Outcomes Report EO.TRA.PCS001.v1. 2021.
8. Lötters FJ, van den Bergh JP, de Vries F, et al. Current and future incidence and costs of osteoporosis-related fractures in the Netherlands: combining claims data with BMD measurements. *Calcif Tissue Int*. 2016;98:235–43.
9. Hossain M, Andrew JG. Is there a difference in perioperative mortality between cemented and uncemented implants in hip fracture surgery? *Injury*. 2012;43:2161–2164.
10. Mundi S, Pindiprolu B, Simunovic N, et al. Similar mortality rates in hip fracture patients over the past 31 years. *Acta Orthop*. 2014;85:54–59.
11. Okike K, Chan PH, Paxton EW. Effect of surgeon and hospital volume on morbidity and mortality after hip fracture. *J Bone Joint Surg Am*. 2017;99:1547–1553.
12. Bentler SE, Liu L, Obrizan M, et al. The aftermath of hip fracture: discharge placement, functional status change, and mortality. *Am J Epidemiol*. 2009;170:1290–1299.
13. Tang VL, Sudore R, Cenzer IS, et al. Rates of recovery to pre-fracture function in older persons with hip fracture: an observational study. *J Gen Intern Med*. 2017;32:153–158.
14. Lu Y, Uppal HS. Hip fractures: relevant anatomy, classification, and biomechanics of fracture and fixation. *Geriatr Orthop Surg Rehabil*. 2019;10:2151459319859139.
15. Swart E, Makhni EC, Macaulay W, et al. Cost-effectiveness analysis of fixation options for intertrochanteric hip fractures. *J Bone Joint Surg Am*. 2014;96:1612–20.
16. Ahn J, Bernstein J. Fractures in brief: intertrochanteric hip fractures. *Clin Orthop Relat Res*. 2010;468:1450–2.
17. Yu J, Zhang C, Li L, et al. Internal fixation treatments for intertrochanteric fracture: a systematic review and meta-analysis of randomized evidence. *Sci Rep*. 2015;5:18195.
18. Sambandam SN, Chandrasekharan J, Mounasamy V, et al. Intertrochanteric fractures: a review of fixation methods. *Eur J Orthop Surg Traumatol*. 2016;26:339–53.
19. Anglen JO, Weinstein JN. Nail or plate fixation of intertrochanteric hip fractures: changing pattern of practice. A review of the American Board of Orthopaedic Surgery Database. *J Bone Joint Surg Am*. 2008;90:700–7.
20. Windolf M, Braunstein V, Dutoit C, et al. Is a helical shaped implant a superior alternative to the Dynamic Hip Screw for unstable femoral neck fractures? A biomechanical investigation. *Clin Biomech (Bristol, Avon)*. 2009;24:59–64.
21. Haidukewych GJ, Israel TA, Berry DJ. Reverse obliquity fractures of the intertrochanteric region of the femur. *J Bone Joint Surg Am*. 2001;83:643–50.
22. Sadowski C, Lübbecke A, Saudan M, et al. Treatment of reverse oblique and transverse intertrochanteric fractures with use of an intramedullary nail or a 95 degrees screw-plate: a prospective, randomized study. *J Bone Joint Surg Am*. 2002;84:372–81.
23. Hoffmann S, Paetzold R, Stephan D, et al. Biomechanical evaluation of interlocking lag screw design in intramedullary nailing of unstable pertrochanteric fractures. *J Orthop Trauma*. 2013;27:483–90.
24. Rudman KE, Aspden RM, Meakin JR. Compression or tension? The stress distribution in the proximal femur. *Biomed Eng Online*. 2006;5:12.
25. Premier Healthcare Database White Paper: Data that informs and performs, March 2 2020. Premier Applied Sciences<sup>®</sup>, Premier Inc. Available at: <https://learn.premierinc.com/white-papers/premier-healthcaredatabase-whitepaper>. Accessed June 8, 2021.
26. Lindenauer PK, Pekow P, Wang K, et al. Perioperative beta-blocker therapy and mortality after major noncardiac surgery. *N Engl J Med*. 2005;353:349–61.
27. Centers for Medicare and Medicaid Services. ICD-10-CM/PCS MS-DRG v37.0. Available at: [https://www.cms.gov/icd10m/version37-fullcode-cms/fullcode\\_cms/P0193.html](https://www.cms.gov/icd10m/version37-fullcode-cms/fullcode_cms/P0193.html). Accessed June 8, 2021.