

Preventing sacral pressure ulcer development in the surgical patient population

Timothy A. Brendle, MS, EN, CNOR, NE-BC



The problem

- It is estimated that 23% of Hospital Acquired Pressure Ulcers (HAPUs) are the result of a patient undergoing a surgical procedure.
- Several factors place surgical patients at increased risk of HAPUs, the most significant being the amount of time a bony prominence, e.g. sacrum, is exposed to pressure. The majority of surgical patients are placed in the supine position, exposing the sacrum to pressure.
- Treatment of pressure ulcers is not only costly, approximately 60,000 deaths are attributed to pressure ulcers each year in the US.

The study

- This trial was conducted at the Ashley River Tower, South Carolina, US.
- 143 patients scheduled for a surgical procedure lasting >3 hours and in a supine position received an ALLEVYN LIFE Sacrum dressing in the sacral area.
- The dressing was implemented within the perioperative area as part of an overall pressure ulcer prevention protocol.

The results

- Of the 143 post-assessments only 1 patient was noted to have Stage 1 (persistent redness).
- Initial evidence shows ALLEVYN LIFE having promise as a pressure ulcer prevention technique when used as part of a pressure ulcer prevention protocol.

Preventing Sacral Pressure Ulcer Development in the Surgical Patient Population



Timothy A. Brendle, MS, RN, CNOR, NE-BC



INTRODUCTION

It is estimated that 23% of Hospital Acquired Pressure Ulcers (HAPU) are the result of the patient undergoing a surgical procedure¹. Numerous factors place the surgical patient at greater risk for developing a HAPU: transfer shear, moisture, ambient room and body temperature, fluid-filled warming/cooling devices, patient's BMI and other comorbidities. However, the most significant factor is the amount of time a bony prominence, such as the sacrum, is exposed to pressure². The majority of surgical patients are positioned in the supine position. This is the rationale for this sacral pressure ulcer improvement project (IP).

LITERATURE REVIEW

As far back as 1959, researchers theorized a causal relationship between surgery and the development of post-surgical pressure ulcers (PU)³. Over the decades, numerous pressure-reducing devices have been used in one form or another to minimize pressure ulcer risks, such as air mattresses, viscoelastic foam, gels and/or a combination of the aforementioned. The devices had varying degrees of success and some devices actually increased PU risks⁴. Not all pressure ulcers can be prevented.

DISCUSSION

The average cost to treat one PU is \$38,000 and even more staggering is approximately 60,000 deaths are attributed to PUs each year in the U.S. alone⁵.


Recent evidence-based practice research utilizing a silicone-foam dressing has shown promise in the prevention of PUs². One academic medical center in the Southeast implemented this IP to evaluate the use of these dressings as a preventive measure at reducing PU incidence in its surgical patient population.

For this IP, all patients at the Ashley River Tower (ART) who were scheduled for a procedure lasting > 3 hours and in the supine position received a silicone hydrocellular foam sacral dressing. A preoperative and postoperative assessment tool was used to determine the clinical effectiveness of preventing erythema/skin breakdown over the sacrum. The tool rated the patients skin from intact thru Stage IV. The assessment also tracked the number of days the dressing was in place and ease of re-applying after daily routine skin assessments.

- ◆ Hyperabsorbent padding layer
- ◆ Pressure distribution
- ◆ Moisture/Shear protection
- ◆ Re-adherence after assessment
- ◆ Adherence average 4.5 days, providing protection post-procedure



Pre-intervention: Percent of OR Patients with Hospital Acquired Pressure Ulcers at both the Main and ART OR's 2013-2014



The silicone hydrocellular foam sacral dressing was implemented within the perioperative area as part of an overall PU prevention program. The silicone hydrocellular foam sacral dressing has numerous PU protective and preventative properties which can help prevent the development of PUs in the surgical patient⁶:

CONCLUSION

During the IP 143 patients received the intervention. Of the 143 post-assessments, only one patient was noted to have a Stage I (persistent redness). Further research, to include random controlled trials, are needed to validate the efficacy of the dressing. However, initial evidence shows promise as a PU prevention technique.

ACKNOWLEDGEMENTS

Heather Craven, MS, RN, CMSRN at the Medical University of South Carolina for assistance with the statistical data. Thanks to the Perioperative staff & Treffe Beaupre at the Ashley River Tower who participated in the IP. The author has no declared affiliation that could be perceived as posing a potential conflict of interest in the publication of this poster.

REFERENCES

1. Connor T, Sledge J, Bryant-Wiersma L, Stamm L, Potter P. Identification of pre-operative and intra-operative variables predictive of pressure ulcer development in patients undergoing urologic surgical procedures. *Urologic Nursing* [serial online], September 2010;30(5):289-305. Available from: CINAHL with Full Text, Ipswich, MA. Accessed August 23, 2014.
2. Brindle T, Wegelin J. Prophylactic dressing application to reduce pressure ulcer formation in cardiac surgery patients. *Wound, Ostomy and Continence Nurses Society*. March/April 2012; 39(2):133-142.
3. Papantonio C, Wallop J, Kolodner K. Sacral ulcers following cardiac surgery: incidence and risks. *Advances in Wound Care*. March 1994; 7(2):24-32.
4. Feuchtinger J, de Gie R, Dassan T, Hafens R. A 4-cm thermoactive viscoelastic foam pad on the operating room table to prevent pressure ulcer during cardiac surgery. *Journal Of Clinical Nursing* [serial online], February 2006;15(2):162-167. Available from: CINAHL with Full Text, Ipswich, MA. Accessed August 24, 2014.
5. Fred C, Ford S, Wagner D, Vanbrackle L. Intraoperatively acquired pressure ulcers and perioperative normothermia: A look at relationships. *AORN*. September 2012; 96(3):251-260.
6. Clarke B. Positive patient outcomes: The use of a new silicone adhesive hydrocellular foam dressing for pressure ulcer prevention and treatment. Smith & Nephew Poster Presented at CAET 2013.