

COBLATION[®] chondroplasty is an effective method to treat partial thickness cartilage lesions

On second-look knee arthroscopy, 88% of previously treated partial-thickness cartilage lesions showed no signs of progression



Study design

- Out of 193 consecutive patients receiving COBLATION treatment from a single surgeon for partial-thickness cartilage defects in the knee, 15 (25 lesions) required repeat arthroscopy for recurrent or new injuries
- Lesion location, size, grade and stability were compared after an average of 10.4 months between procedures



Key results

- Mean lesion size decreased from initial procedure to second-look arthroscopy in all lesions (n = 25; Figure 1a) and tibiofemoral joint lesions (n = 14; Figure 1b)
 - No tibiofemoral lesions showed any signs of progression
- The majority of lesions (88%) showed no signs of progressing at second-look arthroscopy, and 56% improved with complete or partial filling



Figure 1. Mean lesion size (mm²)



Conclusion

Only three of the 25 lesions (12%) in this study demonstrated further deterioration of the cartilage defects after COBLATION treatment. Additionally, partial or complete filling was observed in 56% of the treated lesions.



Considerations

- No marrow stimulating procedures were performed during initial chondroplasty
- Only symptomatic patients warranting repeat arthroscopy were included in this study



Study citation

Voloshin I, Morse KR, Allred CD, Bissell SA, Maloney MD, DeHaven KE. Arthroscopic evaluation of radiofrequency chondroplasty of the knee. *Am J Sports Med.* 2007;35:1702-1707.