

Qentry[®] Integration FAQs for Administrators

1. What is Qentry integration?

Integration occurs between the account's PACS and Qentry (Brainlab's cloud-based file transfer and storage system) and allows you to more securely, quickly, and efficiently send patient images to the Smith & Nephew representatives. Once the images are sent in Qentry the rep can access them in TraumaCad to template and save back to Qentry to communicate with the surgeon or the CareTeam.

2. What are the advantages of integration?

Integration gives you more control, security, and efficiency. Qentry employs advanced encryption and access control technologies to ensure all sensitive medical information is protected. Image routing rules can be configured during initial integration set up and images can be automatically delivered to the Smith & Nephew representative with DICOM push, instead of physically giving the images to the rep for each case.

3. What are my options for pushing images to the representative?

Once you have created a Qentry.com account, there are three ways to push images:

1. Qentry.com: Upload/download files quickly from any computer using Qentry web uploader
2. Qentry Desktop: Access PACS for manual upload/download, upload DICOM CDs and individual images
3. Qentry Gateway: Integrate with PACS for automatic upload and manual download for imaging modalities or other DICOM applications that are able to communicate via DICOM Query/Retrieve (Q/R) or DICOM Push

Refer to the Qentry website "Help" section for more details on each option.

4. With Gateway integration how can I limit the images that the Smith & Nephew rep receives?

During the integration process, you will establish permission so the rep will only receive or access the appropriate images. For additional information go to Qentry.com register for a free account, login, click Help and refer to the [Qentry Gateway Software User Guide](#).

5. Is Qentry HIPAA compliant? Where can I find additional information about Qentry's data security and privacy policy?

Yes, Brainlab is committed to comply with applicable rules and regulations of the Health Insurance Portability and Accountability Act (HIPAA), as amended by the Health Information Technology for Economic and Clinical Health Act (HITECH) through the implementation of stringent privacy and security policies and procedures.

For information regarding Qentry's data security and privacy policy, refer to the following documents in Qentry.com's Help section or refer to [Qentry.com/Security](#).

[Security Statement](#)

[Privacy Policy](#)

6. How much does integration support cost?

Integration support is free. The integration support expenses will be paid by Smith & Nephew.

7. How do I initiate technical support for Gateway integration?

Contact your Smith & Nephew representative for more information.

8. How long does it take to push images into Qentry?

The time required to push images into Qentry will vary depending on the image size and the imaging center's internet bandwidth. For example, 20 studies with 2 images each (total of 1.4GB) on a standard internet connection should take about 28 seconds to upload.

9. What options do I have if we already have TraumaCad®?

There are several options depending on your set up with Brainlab. Please contact your Smith & Nephew representative to set up a call with Brainlab technical support.

10. Who should I contact if I have additional questions?

Please contact VISIONAIRE Logistics for more information at VisDigitalTemp@smith-nephew.com or call VISIONAIRE Support at 1-800-262-3536 (when prompted, select option "2").