

 **smith&nephew**  
**ACCU-PASS<sup>®</sup>**  
Suture Shuttle




Form, Function, and Flexibility



With its single-handed design, ergonomic handle, and thumb wheels, the Smith & Nephew ACCU-PASS<sup>◊</sup> Suture Shuttle provides precise control for suture placement and shuttle deployment, reducing operative time and cost.

Eight different tip configurations and three shuttle loading options provide flexibility in both procedure and technique, while the solid one-piece shaft resists bending and breakage. Sterile, single-patient use gives you the confidence of knowing the instrument is sharp. Each ACCU-PASS Suture Shuttle is designed to accommodate a doubled #0, a doubled #1 or a single #2 monofilament suture.

**Suture passing and shuttling for instability, plications and 'side to side' repairs.**



**Ordering Information**

- 7210426 ACCU-PASS Suture Shuttle, Straight
- 7210427 ACCU-PASS Suture Shuttle, Crescent
- 7210425 ACCU-PASS Suture Shuttle, 45°, Upbend
- 72200419 ACCU-PASS Suture Shuttle, 70°, Upbend
- 7210423 ACCU-PASS Suture Shuttle, 45°, Left
- 7210424 ACCU-PASS Suture Shuttle, 45°, Right
- 72200418 ACCU-PASS Suture Shuttle, Big curve
- 72201076 ACCU-PASS Suture Shuttle, J hook
- 72201361 Sterile single packed #1 monofilament material, box of 10