



Tonsillectomy and adenoidectomy with one device

Versatile surgical instruments for the removal of both  
tonsils and adenoids


 **smith&nephew**  
**COBLATION<sup>®</sup>**  
**EVAC<sup>®</sup> 70 XTRA**  
Adenotonsillectomy Wand

# EVAC<sup>◇</sup> 70 XTRA adenotonsillectomy Wand

Versatile surgical instruments for the removal of both tonsils and adenoids

## Key features

- Triple-wire active electrode configuration removes tissue for both tonsillectomy and adenoidectomy
- Integrated saline and suction ports allow for quick and easy operating room setup
- 6 inch, malleable shaft allows improved access to the choanae during adenoidectomy



## Ordering information

### EVAC 70 XTRA Wand

Reference #	Description
EICA5872-01	EVAC 70 XTRA WAND
EICA5874-01	EVAC 70 XTRA HP WAND
EA2020-01	EVAC BENDING TOOL
EC8000-01	COBLATOR <sup>®</sup> II CONTROLLER

### ArthroCare Corporation

7000 West William Cannon Drive  
Austin, TX 78735  
USA

[www.smith-nephew.com](http://www.smith-nephew.com)

Order Entry: 1-800-343-5717  
Order Entry Fax: 1-888-994-2782

© 2015 Smith & Nephew, Inc.  
™Trademark of Smith & Nephew. Reg. US Pat. & TM Office.