

59th Annual
Orthopaedic Research Society Meeting, 2013
San Antonio, Texas



Oxidized Zirconium vs Cobalt-Chromium TKA: Surface roughness of retrieved femoral components

Marcella E. Elpers, Thomas J. Heyse, Danyal H. Nawabi, Timothy M. Wright, Steven B. Haas
Hospital for Special Surgery, New York, NY

Introduction

- Oxidized zirconium (OxZr) was introduced as an alternative bearing surface to CoCr alloys in TKAs
- Simulator studies have shown that the surface roughness of OxZr does not change during testing, while the surface roughness of CoCr increases¹
- Retrieval studies have shown that PE components paired with OxZr have less surface damage than those paired with CoCr femoral components²
- While in vivo performance of CoCr is well established, little is known about the in vivo performance of OxZr³


Objective

To evaluate the surface characteristics of OxZr and CoCr retrieved TKA femoral components

Methods





Qty – 10 OxZr
Retrieved femoral components



Qty – 10 CoCr
Retrieved femoral components


Implants matched on patient demographics: LOI, Age, Revision Diagnosis and BMI

Surface roughness measured with non-contact white light profilometer
(MicroXAM Optical Profiler, ADE Phaseshift, AZ)




Statistical Analysis: General Estimating Equations (GEE) models were developed to estimate the differences in roughness between the materials.

Results



- Scratching visible on CoCr retrievals, while little is seen on retrieved OxZr components
- The CoCr implants roughened significantly more in vivo than the OxZr components
- The average surface roughness for the retrieved CoCr was 83% greater than the OxZr retrieved components

3D surface profiles of Oxidized Zirconium and Cobalt-Chromium femoral components



* Median values reported for roughness parameters. All values in μm .

Conclusions

- Increased surface roughness with CoCr could impact the *in vivo* adhesive and abrasive PE wear mechanisms
- Understanding the surface roughness of the bearing material can help explain the damage seen *in vivo*
- A smoother surface, such as OxZr, with superior wear properties is an ideal bearing surface for a TKA

References

- 1 White SE. CORR. 1994.
- 2 Heyse, T. Knee. 2011.
- 3 Good, V. JBJS. 2003.

Acknowledgements

Clark and Kirby Foundations, and The Trump Institute for Implant Analysis. This work made use of the Cornell Center for Materials Research Shared Facilities which are supported through the NSF MRSEC program (DMR-1120296).

Smith & Nephew, Inc.
7135 Goodlett Farms Parkway
Cordova, TN 38016
USA

www.smith-nephew.com

Telephone: 1-901-396-2121
Information: 1-800-821-5700
Orders/Inquiries: 1-800-238-7538