







Ordering information

VERSAJET[®] II Hydrosurgery System







Order No	Description
66800039	Console (115V/230V)
66800979	Surgical cart
66800475	Replacement shelf (Retrofit to 50800)

VERSAJET[®] II Exact Handsets for maximum dermal preservation

Order No	Description
 66800041	VERSAJET II Exact disposable handset (45°/14mm)
 66800042	VERSAJET II Exact disposable handset (45°/8mm)
 66800040	VERSAJET II Exact disposable handset (15°/14mm)
	0.11mm Nozzle orifice
	Lower deck height
	Narrower channel

Order No	Description
66800474	Console user manual
66800472	Multi-function footswitch
66800139	Replacement power cord (N. America)

VERSAJET[®] II Plus Handsets for maximum removal of non-viable tissue and contaminants[†]

Order No	Description
 66800044	VERSAJET II Plus disposable handset (45°/14mm)
 66800045	VERSAJET II Plus disposable handset (45°/8mm)
 66800043	VERSAJET II Plus disposable handset (15°/14mm)
	0.10mm Nozzle orifice
	Higher deck height
	Wider channel

[†] The VERSAJET II Plus handpiece is more powerful than the VERSAJET II Exact. VERSAJET II Plus will select, excise and evacuate tissue faster and more aggressively.



Wound Management

Smith & Nephew, Inc.
970 Lake Carillon Drive
Suite 110
St. Petersburg, FL 33716
USA

Customer Care Center
1-800-876-1261
T 727-392-1261
F 727-392-6914

www.smith-nephew.com
www.versajet.info