

# Femoral Heads Neck Length Translator



Size	Europe		US
28mm			-3
	S		+0
	M		+4
	L		+8
	XL	not equivalent to	+12 <sup>1)</sup>
	XXL	not equivalent to	+16 <sup>1)</sup>
32mm			-3
	S		+0
	M		+4
	L		+8
	XL		+12
	XXL		+16
36mm			-3
	S		+0
	M		+4
	L		+8
	XL		+12

<sup>1)</sup> As result of 3.5mm increments on 28mm heads in Europe which cumulates in XL and XXL to a significant difference.

For combinations with components of different manufactures, the surgeon must verify if the combination is allowed as listed in the compatibility matrix (lit. no. 04758) available on [www.smith-nephew.com/compatibilitymatrix](http://www.smith-nephew.com/compatibilitymatrix). Ball heads should only be combined with hip stems of identical taper dimensions.