



Rapid Ablation and Classic Design

 **smith&nephew**
STARVAC[◇]
COBLATION[®] Wand



STARVAC[◇] COBLATION[◇] Wand

The STARVAC 90 COBLATION[◇] Wand with integrated cable is the newest 90° suction wand offered by SMITH & NEPHEW Sports Medicine. The innovative screen design combines aggressive ablation speed with powerful suction. The classic 90° Wand tip profile is optimized for easy access and precise ablation and hemostasis.

Key features

- Innovative single aperture, star screen design enhances suction flow
- Fast, accurate ablation for efficient procedures
- Same classic size, shape, and profile as best-selling SUPER TURBOVAC[◇] 90
- Available with enhanced COAG mode for improved hemostasis when used with the QUANTUM[◇] Controller platform

Suggested indications

- Subacromial decompression
- Bursectomy
- ACL debridement
- Synovectomy


Technical specifications

STARVAC COBLATION Wand

CATALOG #	ASC4251-01
CONTROLLER PLATFORM	QUANTUM OR ATLAS [◇]
TIP DIAMETER SIZE (MAX)	5.5MM
SHAFT SIZE	3.5MM
ANGLE	90°
RECOMMENDED CONTROLLER SETTINGS	7-9
INTENDED USE	SOFT TISSUE ABLATION
ELECTRODES	SCREEN 1, BALL 4



STARVAC 90 COBLATION Wand tip



STARVAC 90 COBLATION Wand screen

Arthrocare Corporation

7000 West William Cannon Drive
Austin, TX 78735
USA

www.smith-nephew.com

Order Entry: 1-800-343-5717
Order Entry Fax: 1-888-994-2782

© 2015 Smith & Nephew, Inc.

[◇]Trademark of Smith & Nephew. Reg. US Pat. & TM Office.

P/N 31632 Rev. C 01/15