

The ease of versatility

 **smith&nephew**  
**LEGION<sup>®</sup>**  
Total Knee System




# Versatility

## Maximize results and achieve peak performance

The LEGION® Total Knee System empowers you to address diverse surgical challenges with one A-P Resection, one system and one company. You can personalize joint care no matter what the situation calls for:


### Bi-cruciate retaining partial knee to total knee arthroplasty

- Quickly convert from the JOURNEY® DEUCE® system




### CR to PS

- Easily add PS Box without having to change A-P Resection
- Supports PCL Management




### PS to PS Augmented

- Quickly adjust for femoral defects




### PS to PS Constrained

- Add a constrained insert without using a stem, if desired



### Primary to Revision



## Femoral options



CR



PS



CK

All femoral components are available in OXINIUM® material or cobalt chrome

---

## Inserts

The LEGION® TKS high flexion inserts provide *up to 155° range of motion*, and are available in CPE and XLPE. Constraint may be dialed in from Cruciate Retaining up to 2° – 3° of varus/valgus constraint.



Cruciate  
Retaining



Deep Dished



CR  
High Flexion



PS  
High Flexion



Posterior  
Stabilized



Constrained\*

## Tibial base plates

Achieve maximum coverage while maintaining proper axial alignment with the Smith & Nephew asymmetrical tibial base plates. Engineered with the clinically proven locking mechanism, these are available in both primary and revision designs.



Primary



Primary



Revision

## Patellas



Oval Resurfacing\*



Round Resurfacing\*



Biconvex\*

Available in 9mm and 13mm thicknesses

\*Not available in XLPE.

Smith & Nephew, Inc.  
7135 Goodlett Farms Pkwy  
Cordova, TN 38016  
USA

[www.smithnephew.com](http://www.smithnephew.com)  
[www.legionpower.com](http://www.legionpower.com)

Telephone: 1-901-396-2121  
Information: 1-800-821-5700  
Orders/inquiries: 1-800-238-7538