



New!

Gentle on tissue
Tough under pressure

Combine with our new CROSSTRAC[®] Hip Guide System for accurate trajectories and first-time access to the hip joint.

 **smith&nephew**
ARTHROGARDE[®]
Hip Access Cannulas



Introducing the new guard for arthroscopy: ARTHROGARDE[◊] Hip Access Cannulas are designed to be atraumatic, featuring smooth ends for optimal safety and high-strength design for enhanced durability.

The high-strength, atraumatic design makes the ARTHROGARDE Hip Access Cannula ideal for hip arthroscopy. The uniquely designed rounded ends provide a smooth surface to protect soft tissue when performing hip arthroscopy. Plus, they're extra strong to prevent bending or wear and tear over time. Combine with our CROSSTRAC[◊] Hip Guide System* for accurate trajectories and first-time access to the hip joint.



Rounded Ends

Rounded ends provide a smooth surface that is atraumatic in design



Connector

High-strength design provides greater durability and stable connections so accurate trajectories can be achieved when using the CROSSTRAC Hip Guide



Slotted Cannula*

Featuring the new atraumatic design, enhanced durability, angled profile and ergonomic handle for easy maneuverability

Ordering Information

ARTHROGARDE Hip Access Cannulas



72201741 4.5 mm ARTHROGARDE Hip Access Cannula GREEN (CROSSTRAC Hub)



72201742 5.0 mm ARTHROGARDE Hip Access Cannula BLUE (CROSSTRAC Hub)



72201743 5.5 mm ARTHROGARDE Hip Access Cannula RED



72201590 ARTHROGARDE Slotted Hip Access Cannula

CROSSTRAC Hip Guide

72201832 CROSSTRAC Hip Guide System

System Includes:

72201831 CROSSTRAC Hip Guide with Bullet

72201741 4.5 mm ARTHROGARDE Hip Access Cannula GREEN (CROSSTRAC Hub)

72201742 5.0 mm ARTHROGARDE Hip Access Cannula BLUE (CROSSTRAC Hub)

72200096 4.5 mm Cannulated Obturator with Delrin Handle

72200097 5.0 mm Cannulated Obturator with Delrin Handle

Sold Separately:

72201811 XL (extra-long) Disposable Arthroscopy Needle