

Stem and Femoral Ball Head Combinations
REDAPT°, SPECTRON°, BASIS° and ECHELON°

✓	Combination is approved by Smith & Nephew
✗	Combination is not approved by Smith & Nephew
□	Combination under consideration

Combinations listed in the matrix may not be approved in all individual markets or geographies. The information contained in this matrix does not supersede the instructions for use (Package insert) in force in the markets where the products are being used. Please refer to your local Smith & Nephew representative to confirm the approval status in your country or region if you have questions about how Smith & Nephew products can be used.

		REDAPT	SPECTRON	BASIS	ECHELON				
		All sizes	All sizes	All sizes	All sizes				
Ball Head Ø	Neck length								
Oxinium	D22	S / +0	M / +4	L / +8	XL / +12	✓	✓	✓	✓
	D26	S / +0	M / +4	L / +8	XL / +12	✓	✓	✓	✓
	D28 to D36	XS / -3				✓	✓	✓	✓
	D28 to D36	S / +0	M / +4			✓	✓	✓	✓
	D28 to D36			L / +8		✓	✓	✓	✓
	D28				XL / +12	✓	✓	✓	✓
	D28					XXL / +16	✓	✓	✓
	D32 and D36				XL / +12	✓	✓	✓	✓
	D32					XXL / +16	✓	✓	✓
	D40 and D44 ³	XS / -4				✓	✓	✓	✓
	D40 and D44 ³		S / +0	M / +4	L / +8		✓	✓	✓
	CoCrMo	D22	S / +0	M / +4	L / +8	XL / +12	✓	✓	✓
D26		S / +0	M / +4	L / +8	XL / +12	✓	✓	✓	✓
D28 to D36		XS / -3				✓	✓	✓	✓
D28 to D36		S / +0	M / +4			✓	✓	✓	✓
D28 to D36				L / +8		✓	✓	✓	✓
D28					XL / +12	✓	✓	✓	✓
D28						XXL / +16	✓	✓	✓
D32 and D36					XL / +12	✓	✓	✓	✓
D32						XXL / +16	✓	✓	✓
D40 and D44 ³		XS / -4				✓	✓	✓	✓
D40 and D44 ³			S / +0	M / +4	L / +8		✓	✓	✓
D22				M	L		✓	✓	✓
D28 and D32		S	M	L		✓	✓	✓	✓
D28 and D32					XL	✓	✓	✓	✓
D28 and D32						XXL	✓	✓	✓
Stainless Steel	D22	S / +0	M / +4	L / +8	XL / +12	✗	✗	✗	✗
	D26	S / +0	M / +4	L / +8	XL / +12	✗	✗	✗	✗
	D28 and D32	XS / -3				✗	✗	✗	✗
	D28 and D32	S / +0	M / +4	L / +8		✗	✗	✗	✗
	D28 and D32				XL / +12	✗	✗	✗	✗
	D28 and D32					XXL / +16	✗	✗	✗
BioloX Forte	D28 to D36	S / +0	M / +4	L / +8		✓	✓	✓	✓
	D28 to D36	S	M	L		✓	✓	✓	✓
BioloX Delta	D32 and D36	S / +0	M / +4	L / +8		✓	✓	✓	✓
	D36				XL / +12	✓	✓	✓	✓
	D40	S / +0	M / +4	L / +8		✓	✓	✓	✓
	D28 to D36	S	M	L		✓	✓	✓	✓
BioloX* Option	D32 and D36				XL	✓	✓	✓	✓
	D28 to D36	S	M	L		✓	✓	✓	✓
BioloX* Option	D28 to D36				XL	✓	✓	✓	✓

³ 40mm and 44mm heads include Ti-sleeves.

This matrix is only applicable for Smith & Nephew Femoral Ball Heads and Stems as well as third-party products (identified by *) covered in this document.

* BioloX Option (Revision THA) Femoral Ball Heads are products of CeramTec GmbH and compatible with Smith & Nephew stems in respect to this matrix. This product is only available for sale in the EU.

This compatibility matrix can only be accessed online via Smith & Nephew's website <http://www.smith-nephew.com/compatibilitymatrix>. It is the responsibility of the user to consult the Smith & Nephew website to ensure the currency of compatibility information.