


Inventory Control Sheet

Patient name:

Date:

To order: 1-800-238-7538



Half Pins

Catalog #	Description	Qty
4.5mm Threaded HA HALF PINS (SST)		
7107-0823	10mm x 150mm	
7107-0824	20mm x 150mm	
7107-0825	30mm x 150mm	
7107-0826	40mm x 150mm	
7107-0827	50mm x 150mm	

Catalog #	Description	Qty
6mm Threaded HA HALF PINS (SST)		
7107-0800	20mm x 150mm	
7107-0801	20mm x 200mm	
7107-0802	25mm x 150mm	
7107-0803	25mm x 200mm	
7107-0804	30mm x 150mm	
7107-0805	30mm x 200mm	
7107-0806	35mm x 150mm	
7107-0807	35mm x 200mm	
7107-0808	40mm x 150mm	
7107-0809	40mm x 200mm	
7107-0810	50mm x 150mm	
7107-0811	50mm x 200mm	
7107-0812	55mm x 150mm	
7107-0813	55mm x 200mm	
7107-0814	60mm x 150mm	
7107-0815	60mm x 200mm	
7107-0816	65mm x 150mm	
7107-0817	65mm x 200mm	
7107-0818	70mm x 150mm	
7107-0819	70mm x 200mm	
7107-0820	80mm x 200mm	
7107-0821	90mm x 200mm	
7107-0828	15mm x 150mm	
7107-0829	15mm x 200mm	
7107-0835	45mm x 150mm	
7107-0836	45mm x 200mm	

Catalog #	Description	Qty
Hex Fix Pins (TI)		
122700	4mm x 20mm	
122701	4mm x 25mm	
122702	4mm x 30mm	
122703	4mm x 35mm	
122704	4mm x 40mm	
122707	5mm x 20mm	
122708	5mm x 25mm	
122709	5mm x 30mm	
122710	5mm x 35mm	
122711	5mm x 40mm	
122712	5mm x 45mm	
122713	5mm x 50mm	
122714	5mm x 55mm	
122715	5mm x 60mm	
7102-1022	6mm x 20mm	
7102-1023	6mm x 25mm	
7102-1024	6mm x 30mm	
7102-1025	6mm x 35mm	
122725	6mm x 40mm	
122727	6mm x 50mm	
122729	6mm x 60mm	

Catalog #	Description	Qty
6mm Threaded HA HALF PIN, Cannulated		
7107-0822	50mm x 150mm Recommended to use 1.6 mm guide wire (7105-1039) or smaller	

Drills

Catalog #	Description	Qty
HA HALF PIN Drills, 4.8 mm Shank (for use with orange handled drill guide 7106-7310)		
7107-0345	3.8mm Stepdrill, short	
7107-0346	4.8mm Drill, short	
7107-0348	3.8mm Stepdrill, long	
7107-0349	4.8mm Drill, long	

Catalog #	Description	Qty
HA HALF PIN Drills, 6mm Shank		
7107-0844	4.8mm, long (use with 6mm pins)	
103046	4.8mm, short (use with 6mm pins)	
103045	3.8mm (use with 4.5mm pins)	
7107-0845	Cannulated 4.8mm (use with 6mm pins). Recommended to use 1.6mm guide wire (7105-1039) or smaller	
7107-0846	Cannulated 3.8mm (use with 4.5mm pins). Recommended to use 1.4mm guide wire or smaller (128042)	

Catalog #	Description	Qty
Hex Fix Drills, 6mm Shank		
103044	2.7mm for 4mm Pins	
103045	3.8mm for 5mm Pins	
103046	4.8mm for 6mm Pins	