



 **smith&nephew**
JET-X[®]
Bar Tibia Kit

Reduction made easy


The JET-X Bar Tibia Kit offers a simple solution for the treatment of tibial fractures.



JET-X[◇] Bar Tibia Kit, 350 mm Bar

The JET-X Bar Tibia Kit is designed for ease of application. The JET-X Bar to Pin Clamp facilitates fracture reduction while providing stability in a low profile and lightweight design.

Design Features



Ball joints angulate 25° allowing an extra degree of freedom to reduce fractures in multiple planes

Lightweight, radiolucent carbon fiber bars

Clamps have independent locking mechanisms, allowing reduction after the frame is attached to the pins

JET-X Bar Tibia Kit, 350 mm Bar

Cat. No. 7106-2705

Includes:

Description	Qty
Bar to Pin Clamp MR Safe	4
350 mm Bar	1
5 mm X 40 mm TiN Half Pin	5
Protective Cap, Black 5 mm	4
Drill Sleeve with Trocar	1
3.5 mm step Drill/AO Conn	1
Pin Driver with 10 mm Hex	1
10 mm Disposable Wrench	1



Sterile Kit for JET-X Bar Tibia Kit