



 smith&nephew

JET-X[®]

Bar Ankle Spanning Kit

Reduction made easy

The JET-X Bar Ankle Spanning Kit offers a simple solution for the treatment of tibial fractures.



JET-X[◇] Bar Ankle Spanning Kit

The JET-X Bar Ankle Spanning Kit is designed for ease of application. The JET-X Bar to Pin Clamp and JET-X Bar to Bar Clamp facilitates fracture reduction while providing stability in a low profile and lightweight design. The ball joint feature provides the ability to easily position the ankle in multiple planes to aid in fracture reduction.

Design Features



Ball joints angulate 25° allowing an extra degree of freedom to reduce fractures in multiple planes

Lightweight, radiolucent carbon fiber bars

Clamps have independent locking mechanisms, allowing reduction after the frame is attached to the pins

Cartridge clamp design allows snap-fit assembly and prevents passive release of bar during reduction.

JET-X Bar Ankle Spanning Kit, Single Clamp

Cat. No. 7106-2702

Includes:

| Description | Qty |
|-----------------------------|-----|
| Bar to Pin Clamp MR Safe | 4 |
| Bar to Bar Clamp MR Safe | 2 |
| 250mm Bar | 1 |
| 300mm Bar | 2 |
| 5mm x 50mm Traction Pin | 2 |
| 5mm x 40mm TiN Half Pin | 3 |
| Protective Cap, Black – 5mm | 7 |
| Drill Sleeve with Trocar | 1 |
| 3.5mm Step Drill/AO Conn. | 1 |
| Pin Driver with 10mm Hex | 1 |
| 10mm Disposable Wrench | 1 |



Sterile Kit for JET-X Bar Ankle Spanning Kit