

Foot & Ankle

CONFIDENTIAL

Discover
Our fami
Innovative
• New tr
• Ancho
• New s
radolu
For more
DYNOMI
2.0 anch
easier a
new des
joint anc
rotating
pre-load
screw-ir
PEEK-OF
Smith &
wrist lip


 **smith&nephew**
RAPTORMITE[◇] 3.0 PK
RAPTORMITE[◇] 3.7 PLA
Suture Anchor with Needles

The mystery of easy deployment has finally been solved



A simple twist—that's all it takes to deploy the RAPTORMITE[◇] Suture Anchor with Needles. With the ideal combination of anchors, sutures and needles for your foot and ankle repairs, the RAPTORMITE Suture Anchor, is big news for small joints.

Gather your own intelligence.



Latest materials and suitable sizes

- 3.0 PEEK-OPTIMA[®] radiolucent, non absorbable anchor
- 3.7 PLLA radiolucent, absorbable anchor

Reliable suture fixation with low knot profile

- Preloaded with two size 0 ULTRABRAID[®] Sutures and needles
- ULTRABRAID Sutures proven to offer superior knot-break strength²

Fast, easy deployment

- Ergonomically designed solid handle
- Twist-and-release* anchor delivery system

**Patent pending*

Ordering Information

Non-absorbable Anchor

- 72201806 RAPTORMITE 3.0 PK Suture Anchor with two size 0 ULTRABRAID Sutures, and needles
- 72202039 2.6 mm Drill and guide for use with RAPTORMITE 3.0 PK Suture Anchor

Absorbable Anchor

- 72201805 RAPTORMITE 3.7 PLLA Suture Anchor with two size 0 ULTRABRAID Sutures, and needles
- 72202038 3.2 mm Drill and guide for use with RAPTORMITE 3.7 PLLA Suture Anchor

1. Data on file at Smith & Nephew 2009.
 2. Based on 2007 benchmark studies comparing ULTRABRAID Sutures to Ethibond[®] and Fiberwire[®] Sutures

For additional information call +1 800 343 5717 or contact your local Smith & Nephew representative.

Endoscopy
 Smith & Nephew, Inc.
 150 Minuteman Road
 Andover, MA 01810 USA

www.smith-nephew.com
 T +1 978 749 1000
 US Customer Service: +1 800 343 5717
 International Customer Service: +1 978 749 1140

©Trademarks of Smith & Nephew. Certain marks Reg US Patent and Trademark Office. All other trademarks acknowledged.
 ©2009 Smith & Nephew, Inc. All rights reserved.
 Printed in USA 05/09 10600460 Rev. A