



Extremities Positioning System

Easy access to tight places.

Positioning for extreme access.

The SPIDER2 Positioner brings measurable traction and unimpeded access to your ankle, wrist and elbow procedures. Hands-free multi-positioning means easy traction and trouble-free positioning and retraction.



Ankle

- Quick and simple setup
- Measurable traction using the TENET Traction Accessory
- Ideal positioning for any angle of knee flexion, including horizontal
- Unobstructed support from underneath operative foot
- Requires counter traction support (FERKEL™ Thigh Holder or TENET Traction Accessory)

Elbow

- Allows unimpeded access to the elbow joint
- Enables easy positioning of elbow in extended or flexed position
- Delivers full control of supination or pronation position with 90° connector
- Facilitates ideal positioning for supine elbow arthroscopy or open elbow procedures

Wrist

- Offers unimpeded access to the surgical wrist
- Delivers measured traction without weights
- Provides exceptional control of both supination and pronation positioning
- Facilitates ideal positioning for pinning and fracture work



Ankle Positioning



Elbow Positioning



Wrist Positioning

Position yourself for the best start.

Distal activation switch can be easily relocated directly to patient. Switch is attached to sterile drape via sticky foam backing – simply remove and adhere to patient and apply a small amount of pressure to activate and move the SPIDER2 positioner.



Step 1



Step 2

Ordering Information

Reference #	Description
WRIST	
72203299	SPIDER2 Limb Positioner
72203300	Switch drape (case of 20)
72203301	SPIDER2 Battery Pack
72200660	Wrist arthroscopy accessory
7210577	TENET Traction Accessory
72200664	Wrist counter traction post
72201550	Small wrist stabilization kit (case of 5)
72201551	Medium wrist stabilization kit (case of 5)
72201552	Large wrist stabilization kit (case of 5)
ANKLE	
72202932	Ankle arthroscopy accessory (2 units)
7210577	TENET Traction Accessory
72202931	Ankle stabilization kit (case of 6)
72201812	FERKEL™ Thigh Holder
ELBOW	
7210572	Elbow connection bar (2 units)
7210749	90° Connector
72200918	SPIDER Cart
Miscellaneous	
72200918	SPIDER Cart
7210553	T-MAX Cart

Endoscopy
Smith & Nephew, Inc.
150 Minuteman Road
Andover, MA 01810 USA

www.smith-nephew.com
T +1 978 749 1000
US Customer Service: +1 800 343 5717
International Customer Service: +1 978 749 1140

Manufactured by TENET Medical Engineering.
Distributed by Smith & Nephew.

©2011 Smith & Nephew. All rights reserved.
Printed in USA. 04/11 10600859 Rev. A



*Data on file from TENET Medical Engineering, file name: Marketing Justification for SPIDER2.