


Take control.
Get better.



Smith+Nephew

CLEARLENS
Tracking Technology

Get better visibility with CLEARLENS Tracking Technology



- Occlusion tolerant markers help to maintain marker visibility¹
- Faster setup and simplified tray management compared to sphere configuration²
- Disposable tracking arrays with preinstalled markers are ready-to use from the peel pack
- Delivered as a sterile packaged set for one knee replacement procedure

Smith & Nephew, Inc.
1450 Brooks Road
Memphis, Tennessee 38116
USA

www.smith-nephew.com

^oTrademark of Smith+Nephew
All Trademarks acknowledged
©2020 Smith & Nephew, Inc.
24588 V1 04/20

 **BRAINLAB** +



References

1. Data on file with Smith+Nephew. Maximum Virtual Shift of Marker Position Due to Occlusion. Form 04-143: Revision 06, Released March 30, 2011.
2. Data on file with Smith+Nephew. Data from one surgeon experience comparing to KNEE2.6.